



# Ultra-sensitive molecular detection of gene fusions from RNA using ASPYPE

Eleanor Gray, Justyna Mordaka, Efthimia Christoforou, Christina Xyrafaki, Kristine von Bargen, Nicola Potts, Barnaby Balmforth

Biofidelity. Ltd, 330 Cambridge Science Park, Cambridge, CB4 0WN, United Kingdom

## Why detect gene fusions from RNA?


- In NSCLC, effective therapies are available to target gene fusions include Alectinib (**ALK**), Selpercatinib (**RET**), Crizotinib (**ROS1**), and Larotrectinib (**NTRK**)<sup>1</sup>
- RNA is an underused analyte, which yields insight into gene expression, and reflects mutations outside of exons
- Single-target gene target assays can result in sample exhaustion due to sequential analysis to cover all targets
- NGS is an important discovery tool when first line screening yields no result ... BUT has long turn-around time, high cost, high sample input requirement, and clinical implications are difficult to interpret

## The ASPYPE assay

- Simple, rapid, low-cost**
- Simple: pre-made mastermixes
- Simple: workflow requires sample input and reagent transfer
- Rapid: **< 1 day to result**
- RNA & DNA analysed simultaneously** in one assay run
- 36 gene fusions** detected
- 22 DNA & 2 RNA PCR wells per patient
- 16 samples per qPCR run

## Equivalence to reference standards

In-house panel & Seraseq® Fusion RNA Mix v4, quantified by dPCR, tested using ASPYPE-lung RNA panel




- Five RNA gene-fusion targets
- Synthetic oligonucleotides or Seracare mix tested
- 0, 3 or 6 copies per sample
- 12 repeats per input level

ASPYPE-lung assay can detect RNA gene fusions across different 3' driver gene families

## Sensitivity - Detection of single-copies

Limiting dilution of sample input to estimate limit of detection



- Four RNA fusion targets
- Synthetic oligonucleotides
- 0, 1, 2, 3, 6 or 9 copies per sample
- 12 repeats per input level

Theoretical number of copies in each reaction	Expected number of positive reactions (out of 12)	ROS1 (CD74-ROS1) positive reactions	ALK (EML4-ALK) positive reaction	RET (TRIM33-RET) positive reactions	NTRK (ETV6-NTRK3) positive reactions
0	0	0	0	0	0
1	7.6 ± 1.7	9	10	9	8
2	10.4 ± 1.2	11	12	12	11
3	11.4 ± 0.8	11	12	12	12
6	12.0 ± 0.2	12	12	12	12
9	12.0 ± 0.03	12	12	12	12

Detection of RNA gene fusion targets by ASPYPE-lung is consistent with single molecule detection limits when compared to expected Poisson distribution for input copy number

## ASPYPE-Lung assay with FFPE Patient Samples: Specific & Sensitive


Specificity - test ASPYPE-lung RNA panel on five non-cancerous lung samples

- Five FFPE lung tissue blocks from patients without diagnosed NSCLC
- Four curls taken per block, and RNA extracted
- 1 ng, 5 ng, and 10 ng RNA tested per curl

100% specificity for ALK, ROS1, RET, NTRK1, NTRK2 and NTRK3 fusions

Sensitivity - test ASPYPE-lung RNA panel on known fusion-positive FFPE tissue at different input levels

- RNA extracted from two fusion-positive patient FFPE lung tissue samples
- Tested using ASPYPE-lung RNA panel (operator blind)
- 1, 2, 5, 10 ng input RNA
- 8 repeats per input sample level



data look suspicious due to bug in code that makes max 127 - needs reanalysing plus thresholds recalculating

Thresholds difficult to pick for this target (based on Lung RNA? Water? Non-NSCLC FFPE?) Once I've got reanalysed data I'll work out what to use

Sample	ROS1 result (reactions)	ALK result (reactions)	RET result (reactions)	NTRK result (reactions)
NTC	Negative (8/8)	Negative (8/8)	Negative (8/8)	Negative (8/8)
Human Lung RNA	Negative (8/8)	Negative (8/8)	Negative (8/8)	Negative (8/8)
ROS1 synthetic target	Positive (8/8)	Negative (8/8)	Negative (8/8)	Negative (8/8)
ALK synthetic target	Negative (8/8)	Positive (8/8)	Negative (8/8)	Negative (8/8)
RET synthetic target	Negative (8/8)	Negative (8/8)	Positive (8/8)	Negative (8/8)
NTRK synthetic target	Negative (8/8)	Negative (8/8)	Negative (8/8)	Positive (8/8)
NSCLC_151: 1 ng, 5 ng, 10 ng	Positive (8/8) at 1 ng, 2 ng, 5 ng	Negative (8/8)	Negative (8/8)	Negative (8/8)
NSCLC_152: 1 ng, 5 ng, 10 ng	Negative (8/8)	Positive (8/8) at 1 ng, 2 ng, 5 ng	Negative (8/8)	Negative (8/8)



One sample was ALK fusion-positive, and one sample was ROS1 fusion-positive by orthogonal testing; ASPYPE-lung RNA results were concordant for both samples

Both samples yielded positive results at all sample input levels, including 1 ng

## ASPYPE assay design - DNA & RNA combined


- ASPYPE-lung panel detects all NCCN-recommended biomarkers, including:

- deletions, insertions and SNVs in EGFR, BRAF, KRAS, ERBB2 (from DNA)<sup>2</sup>
- gene fusions: ALK, RET, ROS1, NTRK1, NTRK2, NTRK3 (from RNA)



## ASPYPE assay output

- Binary result: positive (mutation present) or negative (mutation absent) for panel targets e.g. ALK-positive, EGFR T790M-positive
- Clinical result of mutation present or absent based on fluorescence signal that is higher or lower than a control-based threshold



The raw ASPYPE assay output for a single mutation - interpreted for the user as negative or positive if above or below the separation threshold

## Conclusions

- ASPYPE-lung RNA panel detects 36 gene fusions from RNA
- Panel includes most common ALK, ROS1, RET and NTRK mutations
- Assay workflow takes < 1 day and is run concurrently with DNA sample
- Detection consistent with single molecule detection limits
- 100% sensitivity and specificity from clinical samples across all variants and input quantities

## References

- Tan AC, Tan DSW. Targeted therapies for Lung cancer patients with oncogenic driver molecular alterations. *Journal of Clinical Oncology* 40, no. 6 (February 20, 2022) 611-625. doi: 10.1200/JCO.21.01626.
- Silva AL, Powalowska PK, Stolarek M, Gray ER, Palmer RN, Herman B, Frayling CA, Balmforth BW. Single-copy detection of somatic variants from solid and liquid biopsy. *Scientific Reports* 2021 Mar 16;11(1):6068. doi: 10.1038/s41598-021-85545-3.

